

UNITRON®



BALL BEARING
BOOM STAND

MANUAL

UNPACKING AND COMPONENTS

- Horizontal pillar with cross block
- Base
- Vertical post
- Knobs (3)
- Safety lock collar
- Hex bolt
- Lock lever
- Hex wrench

Warnings: Some components (especially the base) are very heavy, please handle with extra caution.

INSTRUCTIONS

1. Carefully remove all parts from box and check for all package contents as listed above and shown below.



BALL BEARING BOOM STAND

2. Place the base on the edge of a table as shown below. Use the hex bolt to attach the vertical post to the base, and tighten using the hex wrench provided.



3. Place the base on the table and unscrew the cap on the vertical post.



BALL BEARING BOOM STAND

4. Insert the safety lock collar onto the vertical post and tighten/lock it with one of the knobs.



5. Slide the horizontal pillar with cross block into the vertical post.



BALL BEARING BOOM STAND

- Adjust to the desired position and tighten the cross block on both sides using the other 2 knobs.

Put the cap back into the vertical post.



- Screw the lock lever on the end of the horizontal pillar. This lock lever is to secure the e-arm or focus mount.



BALL BEARING BOOM STAND

8. Assembly is complete.




Note: This ball bearing boom stand is compatible with all UNITRON® stereo microscopes and accessories.

E-ARM or FOCUS MOUNT SETUP

1. Insert the e-arm or focus mount into the end of the horizontal bar and tighten using the lock lever.



2. Adjust the e-arm or focus mount to desired level or position.




BALL BEARING BOOM STAND

3. E-arm or focus mount can be inserted from the TOP, or BOTTOM or HORIZONTALLY of the horizontal bar. Set e-arm or focus mount up at various angles and tighten the lock lever for oblique viewing.



TOP MOUNTED



BOTTOM MOUNTED



HORIZONTALLY MOUNTED

v020725