

Outstanding optical
performance and
160mm of height
adjustment for
enhanced ergonomics



Ergo Trinocular Head

for Z6/Z8/Z10 Stereo Microscopes

UNITRON®

The Ergo Trinoc Head . . .

- Provides a market-leading 160mm of binocular eyetube height adjustment
- Allows attachment of digital camera to the viewing head (requires C-mount adapters, not included)
- Mounts easily onto the UNITRON Z6, Z8 and Z10 mag changers
- Is compatible with the Nikon SMZ800N, SMZ1270 and SMZ1270i models
- Choose from 50:50 (CAT# 111-11-05) or 0:100 (CAT# 111-11-06) light split to camera port



Ergo Trinocular Head

(CAT# 111-11-05)

*Shown with optional eyepieces
and 1.0x C-mount adapter*



0.5x C-mount adapter
(CAT# 111-25-05A)



1.0x C-mount adapter
(CAT# 111-25-01A)



0.7x C-mount adapter
(CAT# 111-25-07)



**Ergo trinocular head shown with a
Z10, eyepieces, pole stand,
Excelis™ HD camera,
and 0.5x C-mount adapter**
(All sold separately)

UNITRON®

73 Mall Drive • Commack, NY 11725 • 631-543-2000 (P)
Email: info@unitronusa.com • www.unitronusa.com

