

Bulb Replacement Guide

Some of our microscope models are designed with specific stands, while others can be configured on any number of our stands.

Be sure to select the bulb for your specific microscope or stand.

If you have any questions, contact your authorized UNITRON dealer, or call us at 631-543-2000.

Model/Catalog #	Bulb Catalog #	Bulb Specifications
EXAMET-4 Reflected & Transmitted	16049	12v 50w halogen transmitted (bottom)
FSB	15887	Fluorescent ring bulb (4 wire lead type)
	15887.1	Fluorescent ring bulb (4 pin type)
MEC2 Microscope	149-30-30	6v 30w halogen bulb
NEOMET	15762	8v 15w EL-38 Long Life (T) bulb
RMM2 Microscope	15435	6v 20w halogen bulb
VERSAMET 2	16049	12v 50w halogen bulb
VERSAMET 3	16048	12v 100w halogen bulb
ZSB	15887	Fluorescent ring bulb (4 wire lead type)
	15887.1	Fluorescent ring bulb (4 pin type)
Unilite Ring Light-70mm/#12884	12885	10w fluorescent bulb
Unilite Ring Light-60.8mm/#12886	12887	10w fluorescent bulb
Diascopic Transfer Stand	141-30-T	12v 15w halogen
	6053A	OR 21v 150w fiber optic